




# Lumis 150 ST with Oxygen and non-invasive circuit set up



**Non vented mask**  
**Note: Good mask fit is important to minimize risk of cross contamination**  
**Select mask setting: Nasal or Full Face**

**Anti-viral filter or anti-viral HMEF**

**Max pressure: 25cm H<sub>2</sub>O**  
**Max backup breath rate: 50bpm**  
**Note: no alarms, only Leak Alert**  
 iVAPS mode is **not** recommended for COVID-19

|                                     |   |
|-------------------------------------|---|
| 28125 or 28226                      | Lumis 150 ST or Lumis 150 VPAP ST-A   |
| 26960                               | Oxygen port connector   |
| 21948                               | Standard tube single  |
| 24988                               | Passive exhalation flow   |
| 24966<br>24967<br>24969             | <b>Anti-viral filter (50 pack)</b><br><b>Anti-viral HMEF (Heat Moisture Exchange Filter) for adults (50 pack)</b><br><b>HMEF for paediatric (50 pack)</b><br>*Any anti-viral filter can be used with this circuit |
| Supplied by hospital                | <b>Oxygen tube</b>  |
| S - 60971<br>M - 60972<br>L - 60973 | <b>Non-Vented AcuCare™ Full Face Mask F1-0</b>  |


|       |  |
|-------|--|
| 37302 | <b>Oximetry Module</b>   |
| 22374 | <b>ResMed Xpod 3012 LP oximeter</b>                                |
| 22371 | <b>Nonin Xpod clip</b>   |
| 22338 | <b>NONIN 7000A Single Use (Disposable) Oximetry Sensors - 24pk</b> |

**Note:**

- Patients receiving therapy via a non-invasive circuit should be cared for in airborne isolation rooms and staff should wear full PPE (including N95/P2 masks) while in the patient's room.
- Please ensure the "AB Filter" option is set to "YES" when set up is complete. If using a F&P humidifier, set 'external Humidifier' to 'YES'

**Oxygen note:**

- Lumis can take supplemental oxygen up to 15L/min in all modes except iVAPS (4L/min).
- SpO<sub>2</sub> monitoring available on device
- SpO<sub>2</sub> monitoring on device to be used with Oximetry

 ResMed Ltd  
 1 Elizabeth Macarthur Drive  
 Bella Vista NSW 2153 Australia  
 +61(02) 8884 1000 or 1800 658 189