

Awaken your best.

There's a full face mask
for every sleeper



Drift to sleep wearing our most comfortable gear

Let's find a mask that fits your needs

At ResMed, we believe a mask should seal securely and comfortably – no matter how you sleep or breathe. That's why each of our masks is designed with different features to meet a variety of needs. After all, every person and every face is unique. With the number of options we have available, we're confident you'll discover one that feels tailor-made.



FULL
FACE

What is a full face mask?

A full face mask typically has a cushion that covers both your mouth and nose to create a seal over both airways.

If you sleep and breathe with your mouth open or experience nasal congestion, you'll be recommended a full face mask. They are also more suitable for therapy at higher pressures as they are more tolerable than other mask types.



MONEY BACK GUARANTEE

We're so sure you'll love our AirFit™ full face masks that we are offering a 30-day money-back guarantee if you're not completely satisfied.*

*Terms and conditions apply. See in store for details.

Rest assured, it's easy to get the right mask

You're sorted with our mask categories

Our full face mask range will continue to grow and get better, and with our new mask categories, finding the right mask has never been easier. These categories make it easy for you to understand the key benefits of each mask, so you can take one home with confidence. Each category offers something a little different. Have a look to see which one suits you best.

Minimalist category

Small. Light. Simple.

Freedom category

Sleep and move comfortably.

Ultra Soft category

Patented memory foam offers unique comfort.

Universal Fit category

Classic designs that fit nearly every face.

Speak to your ResMed Authorised Dealer to find the right mask that will help you adjust to therapy, and receive optimum results night after night.

AirFit™ F30i

Full face mask



Freedom

AirFit F30i is an UltraCompact full face mask that helps keep air tubing up, out of sight and out of the way. It's designed to let you move at night and sleep in many positions while achieving a reliable seal.¹ It's also ideal if you're looking for visual freedom in bed.

AirFit F30i is available in S/SML, S/STD, M/STD, W/STD and M/LGE

- 1 **Flexible conduit frame** allows therapy air to flow through the elbow and into the frame and cushion.
- 2 **UltraCompact cushion** sits under the nose, eliminating red marks and discomfort on your nasal bridge.
- 3 **Magnetic clips** guide the headgear to the frame quickly and easily.²
- 4 **Optional soft wraps** can help add softness and comfort at the cheeks.




1 ResMed onsite fitting study of 73 current full face CPAP mask users and 14 sleep professionals from Nov. 5-14, 2018 in multiple US locations who trialed AirFit F30i and a competitor full face mask.

2 Use of masks with magnetic components is contraindicated in patients with the following pre-existing conditions: a metallic hemostatic clip implanted in the head to repair an aneurysm and/or metallic splinters in one or both eyes.

AirFit™ F30

Full face mask




Minimalist

AirFit F30 is a small, light and UltraCompact full face mask. It sits discreetly under your nose, not over it, so you can wear glasses, read or watch TV in bed. This also means there's no risk of waking up with discomfort or red marks on the bridge of your nose.

AirFit F30 is available in Small and Medium sizes.
Compatible with AirMini.

- 1 **QuietAir™ diffuser vent** minimises noise and gently disperses exhaled air.
- 2 **UltraCompact cushion** sits under the nose, eliminating red marks and discomfort on your nasal bridge.
- 3 **Magnetic clips** guide the headgear to the frame quickly and easily.
- 4 **Quick-release elbow** makes it simple to disconnect from the mask tubing without removing the mask.





AirFit™ F20

Full face mask

Universal Fit

Our AirFit F20 full face masks combine the tried and tested format many users know and trust with our best technologies. With a fit range of 96%^{1-3*}, these masks comfortably fit almost all facial shapes and sizes, so you can rely on these masks for performance.

AirFit F20 is available in Small, Medium and Large sizes. AirFit F20 for Her is available in Small and Medium sizes. Compatible with AirMini.



- 1 **InfinitySeal™ silicon cushion** is interchangeable with AirTouch™ F20 UltraSoft™ memory foam cushion.
- 2 **QuietAir™ diffuser vent** minimises noise and gently disperses exhaled air.
- 3 **Magnetic clips** guide the headgear to the frame quickly and easily.
- 4 **Quick-release elbow** makes it simple to disconnect from the mask tubing without removing the mask.

*References on back page.



AirTouch™ F20

Full face mask

Ultra Soft

AirTouch F20 features our UltraSoft memory foam*, which is designed to provide a light and breathable seal. Because the foam is softer than silicone, it may help reduce the risk of red marks and skin irritation.

AirTouch F20 also shares the same frame and headgear as AirFit F20, as well as the QuietAir diffuser vent, magnetic clips and quick-release elbow.

AirTouch F20 is available in Small, Medium and Large sizes.

Compatible with AirMini.



14
DAY

REPLACEMENT GUARANTEE

If you are not completely satisfied with the UltraSoft memory foam cushion, you may exchange the cushions that are a part of the AirTouch F20 Starter Kit for one AirFit F20 InfinitySeal™ silicone cushion within 14 days of purchase.**

*Recommended to be replaced monthly.

** T&C's apply. See in store for details.



BECOME A VIP TODAY...

Join sleepvantage and be treated to a range of exclusive member benefits.*

- ✓ **Members-only warranty****
Additional sleepvantage membership warranty on ResMed products.
- ✓ **Feel supported**
Access to our friendly product support line.
- ✓ **Stay informed**
Access to a range of educational articles and videos.
- ✓ **Get special offers**
Be part of exciting promotions and special offers.

For more information, contact sleepvantage today



sleepvantage.com.au



1300 305 705 (AUS) or 0800 737 633 (NZ)



info@sleepvantage.com.au

*sleepvantage membership is available to residents of Aus and NZ.

** Terms and conditions apply. See website for more information.

AirFit F20 References

1 ResMed AirFit F20 internal Australian fitting study of 27 existing CPAP patients, conducted between March-April 2016.

2 ResMed AirFit F20 internal USA fitting study of 34 existing CPAP patients, conducted April 2016.

3 ResMed AirFit F20 internal EU and APAC fitting study of 90 existing CPAP patients, conducted June 2016.



ResMed Pty Ltd

1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

AirFit, AirMini, AirTouch, InfinitySeal, QuietAir, sleepvantage and UltraSoft are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see ResMed.com/ip. Specifications may change without notice.

© 2022 ResMed Pty Ltd. 10111834r2 2022-01