





ResMed sleep therapy device comparison

ResMed

	 AirSense™ 10 Elite	 AirSense™ 11 Elite	 AirSense™ 10 AutoSet™	 AirSense™ 10 AutoSet for Her	 AirSense™ 11 AutoSet
Operating pressure range (cm H ₂ O)	4–20	4–20	4–20	4–20	4–20
Device default pressure (cm H ₂ O)	Start pressure: 4 Set pressure: 10	Start pressure: 4 Set pressure: 10	Start pressure: 4 Min pressure: 4 Max pressure: 20	Start pressure: 4 Min pressure: 4 Max pressure: 20	Start pressure: 5 Min pressure: 5 Max pressure: 20
Therapy modes					
CPAP	✓	✓	✓	✓	✓
AutoSet			✓	✓	✓
AutoSet for Her				✓	✓
Product features					
CSA detection	✓	✓	✓	✓	✓
RERA detection		✓		✓	✓
CSR detection	✓	✓	✓	✓	✓
AutoRamp™ with sleep onset detection	✓	✓	✓	✓	✓
Easy-Breathe™ waveform with EPR™	✓	✓	✓	✓	✓
SmartStart™/Stop feature	✓	(Independent)	✓	✓	(Independent)
★ Care Check-In		✓			✓
★ Personal Therapy Assistant*		✓			✓
★ Bluetooth® connectivity		✓			✓
★ Test Drive		✓			✓
★ OTA: over-the-air software updates*		✓			✓
AutoSet Response comfort setting			✓	✓	✓
Optional oximetry	✓	(Bluetooth)	✓	✓	(Bluetooth)
Supplemental oxygen use	Up to 4 L/min	Up to 15 L/min	Up to 4 L/min	Up to 4 L/min	Up to 15 L/min
Humidification					
Standard fully integrated	✓	✓	✓	✓	✓
Climate Control with ClimateLineAir™†	✓	✓	✓	✓	✓
Data					
Built-in cellular technology*	✓	✓	✓	✓	✓
SD data card option	✓	✓	✓	✓	✓
myAir™	✓	✓	✓	✓	✓
AirView™ with remote assist*	✓	✓	✓	✓	✓
Data storage					
Detailed data on SD card (30 sessions)	✓	✓	✓	✓	✓
High-resolution flow on SD card (30 sessions)	✓	✓	✓	✓	✓
Compliance and summary data on SD card/device (365 sessions)	✓	✓	✓	✓	✓

ALWAYS READ THE LABEL AND FOLLOW THE INSTRUCTIONS FOR USE. You should speak to your doctor about your symptoms and whether a CPAP device is suitable for you.

* Cellular communication depends on network availability. Some AirSense 10 devices may not have cellular communication.

† ClimateLineAir 11 is only compatible with AirSense 11 devices.

AirSense, AirView, AutoRamp, AutoSet, ClimateLineAir, Easy-Breathe, EPR, myAir, and SmartStart are trademarks and/or registered trademarks of the ResMed family of companies.

For patent information, see ResMed.com/ip. Specifications may change without notice. © 2023 ResMed Ltd. 10112565 2023 - 07

★ = **NEW FEATURE**