



size

A unique business opportunity

AirMini provides a unique business opportunity. Offered as a primary system for the home or an additional system for travel, its versatility means you can target a wider range of new and existing patients. If you're looking to help patients better meet their lifestyle needs and comfort expectations while growing your business, AirMini presents an attractive solution.



- 1 Based on volume comparisons done by ResMed internal testing prior to the time this piece was printed.
- 2 In a typical overnight temperature of 16-22 degrees C and a dry room 20-30% RH, the HumidX provides humid air 70-80% RH to breathe.
- 3 When compared to a best in class CPAP machine with a heated humidifier on a mid-range setting and unheated tubing (e.g. Aironse 10, humidifier setting 4, SlimLine tube).
- 4 Currently available for selected full face, nasal and nasal pillows mask users.

ResMed Ltd
 1 Elizabeth Macarthur Drive
 Bella Vista NSW 2153 Australia

AirFit, AirMini and HumidX are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see ResMed.com/ip. See ResMed.com for other ResMed locations worldwide. © 2020 ResMed Ltd. Specifications may change without notice. 10110024/3 2020-12

Rx Only



AirMini™



Rewrite your bedtime stories

ResMed.com/AirMini



Actual

What sets AirMini™ apart?

In addition to being the smallest portable CPAP on the market¹, there are a number of reasons why AirMini isn't your ordinary CPAP device.



ActiveAir™

This unique technology regulates airflow to keep the device stable and comfortable as pressure increases.



HumidX™

Our HumidX waterless humidification system allows exhaled air to pass through before it reaches the mask's vent. In typical overnight conditions² this helps ensure minimal heat and moisture loss and allows HumidX to deliver a similar level of humidification as a traditional bedside device.³⁻⁴



AirMini App

ResMed's AirMini app lets you set up your patient's AirMini via Bluetooth using yours or your patient's smart device. Once connected to the device, easy to follow prompts allow for quick and easy therapy setup.



A complete system

The AirMini device, optimised masks, AirMini App and its optional accessories have all been designed to work together to provide a smooth therapy experience for your patients.

The revolution you've been waiting for

Never before has a CPAP device so small been able to deliver on waterless humidification or therapy performance in the same way. Leveraging ResMed's most proven algorithms and revolutionary new technologies, AirMini is ideal for patients looking for high quality therapy in a pocket-sized design that can be used at home or while they're away.



A new contemporary design

AirMini is sleep therapy as it should be – sleek, modern, and compact with smart technology. Whether your patients want a device that sits discreetly in the bedroom or one that packs away neatly for travel, AirMini meets their lifestyle needs with ease. More than comfort and performance, AirMini gives patients the confidence to take their therapy anywhere.

Optimised masks for AirMini



AirMini has been designed to work with our latest AirFit™ and AirTouch™ series masks. Patients can choose from a range of full face, nasal and nasal pillows options – each designed to deliver optimal sleep therapy that is comfortable wherever they sleep.

Compatible masks include:



AirFit P10 for AirMini



AirFit N20*



AirFit N20 Classic



AirTouch N20*



AirFit N30 for AirMini



AirFit F20*



AirTouch F20



AirFit F30

*for Her version also available