



## The right mask for you

At ResMed we believe that mask performance is non-negotiable. A mask should seal firmly but comfortably, no matter how you sleep, and it should not be a struggle to find a model that suits you. That's why ResMed has created a wide range of masks that have an array of features designed to meet different needs. After all, every face is different.

The masks featured in this brochure are examples of ResMed's full face range.



Full face mask

### What is a full face mask?

A full face mask typically has a cushion that covers both your mouth and nose to create a seal over both airways. If you breathe with your mouth open or experience nasal congestion, you may need a full face mask. They are also more suitable if you are prescribed higher pressures as they are more tolerable than other types of masks.

Speak to your ResMed Authorised Dealer about finding the right mask to help you adjust to therapy, and receive optimum results night after night.



## Become a VIP today...

...and take advantage of a range of exclusive member benefits.\*



### Extended warranty\*\*

50% warranty extension on ResMed products.

### Feel supported

Access to our friendly product support line.

### Stay informed

Access to a range of educational articles and videos.

### Get special offers

Be part of exciting promotions and special offers.

### For more information, contact sleepvantage today:

🌐 [sleepvantage.com.au](http://sleepvantage.com.au)


☎ 1300 305 705 (Aus) or 0800 737 633 (NZ)

✉ [info@sleepvantage.com.au](mailto:info@sleepvantage.com.au)

\*sleepvantage membership is available to residents of Aus and NZ.  
\*\*Terms and conditions apply. See website for more information.

**ResMed Ltd**  
1 Elizabeth Macarthur Drive  
Bella Vista NSW 2153 Australia  
+61 (2) 8884 1000 or 1800 658 189

AirFit, AirTouch, InfinitySeal, QuickFit, UltraSoft, QuietAir and sleepvantage are trademarks of ResMed Ltd. Specifications may change without notice. © 2019 ResMed Ltd. 10111225 2019-03



# ResMed's full face range of masks



## AirTouch™ F20

Full face mask

Our softest mask ever


AirTouch F20 is the first mask to use UltraSoft™ memory foam\*, crafted for a perfect blend of comfort and performance. It's designed to adapt to the unique contours of each individual face, creating a light, breathable seal that stays in place throughout the night.



If you are not completely satisfied with the UltraSoft memory foam cushion, you may exchange the cushions that are a part of the AirTouch F20 Starter Kit for one AirFit F20 InfinitySeal™ silicone cushion within 14 days of purchase.

AirTouch F20 is available in Small, Medium and Large sizes. AirTouch F20 uses the same frame as the AirFit F20.

\*To ensure optimum comfort and seal, we recommend replenishing the AirTouch F20 cushion every month.



**InfinitySeal™**  
silicon cushion interchangeable with AirTouch F20 UltraSoft memory foam cushion.

**QuietAir™**  
diffuser vent minimises noise and gently dispenses exhaled air.

**Magnetic headgear clips**  
guide the headgear to the frame quickly and easily.

**Quick-release elbow**  
simple to disconnect from the mask tubing without removing the mask.

## AirFit™ F20


Full face mask

Full face benefits in a compact size

AirFit F20 and AirFit F20 for Her full face masks are designed to comfortably fit all facial shapes and sizes, and will appeal to those looking for a less confined mask that delivers a clear line of sight.



AirFit F20 is available in Small, Medium and Large sizes. AirFit F20 for her also available in Small and Medium sizes.



**QuietAir™**  
diffuser vent minimises noise and gently dispenses exhaled air.

**Magnetic headgear clips**  
guide the headgear to the frame quickly and easily.

**Quick-release elbow**  
simple to disconnect from the mask tubing without removing the mask.

## AirFit™ F30

Full face mask

Minimal coverage, maximum performance

AirFit F30 is a full face mask that sits gently under your nose, not over it. With nothing over your nose, red marks or discomfort on the bridge of your nose can be a thing of the past.



AirFit F30 is available in Small and Medium sizes.



We're so sure that you will love our AirFit full face masks, that we are offering a 14 day money back guarantee if you are not completely satisfied.\*

\*Terms and conditions apply