

Live life in high Res

There's a full face
mask for every
sleeper



Drift to sleep wearing our most comfortable gear

Let's find a mask that fits your needs

At ResMed, we believe a mask should seal securely and comfortably – no matter how you sleep or breathe. That's why each of our masks is designed with different features to meet a variety of needs. After all, every person and every face is unique. With the number of options we have available, we're confident you'll discover one that feels tailor-made.




FULL
FACE

What is a full face mask?

A full face mask typically has a cushion that covers both your mouth and nose to create a seal over both airways.

If you sleep and breathe with your mouth open or experience nasal congestion, you'll be recommended a full face mask. They are also more suitable for therapy at higher pressures as they are more tolerable than other mask types.



MONEY BACK GUARANTEE

We're so sure you'll love our AirFit™ full face masks that we are offering a 30-day money-back guarantee if you're not completely satisfied.*

*Terms and conditions apply. See in store for details.

Rest assured, it's easy to get the right mask

You're sorted with our mask categories

Our full face mask range will continue to grow and get better, and with our mask categories, finding the right mask has never been easier. These categories make it easy for you to understand the key benefits of each mask, so you can take one home with confidence. Each category offers something a little different. Have a look to see which one suits you best.

Versatile Fit Masks

Designs that fit nearly every face.

Minimalist Masks

Small. Light. Simple.


Freedom Masks

Sleep and move comfortably.

Speak to your ResMed Authorised Dealer to find the right mask that will help you adjust to therapy, and receive optimum results night after night.


AirFit™ F30i

Full face mask



DID YOU KNOW?

This mask is also available without magnetic clips.



Freedom

AirFit F30i is an UltraCompact full face mask that helps keep air tubing up, out of sight and out of the way. It's designed to let you move at night and sleep in many positions while achieving a reliable seal.¹ It's also ideal if you're looking for visual freedom in bed.

AirFit F30i is available in S/SML, S/STD, M/STD, W/STD and M/LGE. Non-magnetic option available.

- 1 Flexible conduit frame allows therapy air to flow through the elbow and into the frame and cushion.
- 2 UltraCompact cushion sits under the nose, eliminating red marks and discomfort on your nasal bridge.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.^{2,3}
- 4 Optional soft wraps can help add softness and comfort at the cheeks.




1 ResMed onsite fitting study of 73 current full face CPAP mask users and 14 sleep professionals from Nov. 5-14, 2018 in multiple US locations who trialled AirFit F30i and a competitor full face mask.

2 The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.

3 Non-magnetic option available.

AirFit™ F40



Minimalist


AirFit F40 is ResMed's smallest ultra-compact full face mask with a flexible design that allows you to sleep in your favourite positions.

AirFit F40 is available in SW/STD, M/STD and L/STD sizes.

- 1 Ultra-compact AdaptiSeal™ cushion made of a soft, total silicone rests securely under the nose and adapts to your face.
- 2 Flexible curved mask frame sits below the eyes, designed for comfort and to reduce facial marks.
- 3 QuietAir™ diffused venting gently disperses exhaled air for a more peaceful bedroom environment.
- 4 Quick-release elbow with short tube makes it easy to connect and disconnect to therapy while also reducing tube drag.
- 5 Magnetic clips easily snap together for fast fitting.¹



¹ The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.



AirFit™ F20

Full face mask

DID YOU KNOW?

This mask is also available without magnetic clips.

Versatile Fit

Our AirFit F20 full face masks combine the tried and tested format many users know and trust with our best technologies. With a fit range of 96%^{1-3*}, these masks comfortably fit almost all facial shapes and sizes, so you can rely on these masks for performance.

AirFit F20 is available in Small, Medium and Large sizes. AirFit F20 for Her is available in Small and Medium sizes. Compatible with AirMini. Non-magnetic option available.



- 1 InfinitySeal™ silicon cushion is interchangeable with AirTouch™ F20 UltraSoft™ memory foam cushion.
- 2 QuietAir™ diffuser vent minimises noise and gently disperses exhaled air.
- 3 Magnetic clips guide the headgear to the frame quickly and easily.^{4,5}
- 4 Quick-release elbow makes it simple to disconnect from the mask tubing without removing the mask.

*References on back page.



AirTouch™ F20

Full face mask

Versatile Fit

AirTouch F20 features our UltraSoft memory foam*, which is designed to provide a light and breathable seal. Because the foam is softer than silicone, it may help reduce the risk of red marks and skin irritation.

AirTouch F20 also shares the same frame and headgear as AirFit F20, as well as the QuietAir diffuser vent, magnetic clips⁴ and quick-release elbow.

AirTouch F20 is available in Small, Medium and Large sizes.

Compatible with AirMini.



14
DAY

REPLACEMENT GUARANTEE

If you are not completely satisfied with the UltraSoft memory foam cushion, you may exchange the cushions that are a part of the AirTouch F20 Starter Kit for one AirFit F20 InfinitySeal™ silicone cushion within 14 days of purchase.**

*Recommended to be replaced monthly.

** T&C's apply. See in store for details.



BECOME A VIP TODAY...

Join sleepvantage and be treated to a range of exclusive member benefits.*

- ✓ **Members-only warranty****
Additional sleepvantage membership warranty on ResMed products.
- ✓ **Feel supported**
Access to our friendly product support line.
- ✓ **Stay informed**
Access to a range of educational articles and videos.
- ✓ **Get special offers**
Be part of exciting promotions and special offers.

For more information, contact sleepvantage today



sleepvantage.com.au



1300 305 705 (AUS) or 0800 737 633 (NZ)



info@sleepvantage.com.au

*sleepvantage membership is available to residents of Australia and New Zealand.

**Terms and conditions apply. See website for more information.

AirFit F20 & AirTouch F20 References

1 ResMed AirFit F20 internal Australian fitting study of 27 existing CPAP patients, conducted between March-April 2016.

2 ResMed AirFit F20 internal USA fitting study of 34 existing CPAP patients, conducted April 2016.

3 ResMed AirFit F20 internal EU and APAC fitting study of 90 existing CPAP patients, conducted June 2016.

4 The mask contains magnets that may interfere with certain implants or medical devices. Please refer to the User Guide for complete labelling information, including magnet contraindications and warnings.

5 Non-magnetic option available

ALWAYS READ THE LABEL AND FOLLOW THE INSTRUCTIONS FOR USE. You should speak to your doctor about your symptoms and whether a CPAP device is suitable for you.



ResMed Pty Ltd

1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

AirFit, AirMini, AirTouch, InfinitySeal, QuietAir, sleepvantage and AdaptiSeal are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see ResMed.com/ip. Specifications may change without notice.

© 2024 ResMed Pty Ltd. 10111834r4 2024-07